

# Hall Pyke

---



**Membrane Filtration**

**AdvanLife<sup>®</sup>**


**Filter Cartridges**

**The SMART choice for filtration**

# AdvanLife® Filter Cartridges

Highly Asymmetric PES Membrane · Sterilizing-Grade

**AdvanLife®** Filter Cartridges are constructed of a single-layer asymmetric hydrophilic PES membrane. Characteristics include excellent throughput and high dirt hold capacity and durability. The high flow rates, when compared to other sterilizing-grade filter cartridges, reduce filtration costs.



## Features and Benefits

- Highly asymmetric PES membrane provides high dirt holding capacity and longer service life
- Each filter is individually Integrity Tested prior to leaving the factory
- Available in absolute ratings from 0.1µm to 1.2µm for precise bacteria and particle removal
- Complies with Food Contact Regulations: FDA 21CFR177-182 and 1935/2004 EC

## Materials of Construction

Filter Media	PES
Cage/Support	Polypropylene
Core/End Caps	Polypropylene

## Microbiological Stability

Longer life with highly asymmetric membrane






## Operating Conditions

<b>Max. Temperature</b>	80°C
<b>Max. Differential Pressure</b>	4.0 bar / 21°C (forward) 2.4 bar / 70°C (forward)
<b>Bubble Point</b>	BP: >0.34 MPa (water), 0.22 µm BP: > 0.22 MPa(water), 0.45 µm
<b>Diffusion Flow</b>	DF: < 30 ml/min/10" @0.25Mpa(water), 0.22µm DF: < 28 ml/min/10" @0.16Mpa(water), 0.45 µm
<b>Steam Sterilization</b>	121°C/30 min @Max. Differential Pressure for 0.3 bar
<b>Hot Water Sterilization</b>	85°C/30 min @Max. Differential Pressure for 2 bar
<b>Cleaning Solution</b>	2% NaOH Solution @ ≤65°C
<b>Effective Filtration Area</b>	0.58m <sup>2</sup> / Φ69-10 inch

## Flow Rate Characteristics



## Reliable Microbiological Control

The primary purpose of a membrane filter cartridge in beverage processing is to effectively control spoilage microorganisms.

Typical Log Reduction Value (LRV)			
	<b>B.Diminuta</b>	<b>Lactobacillus Brevis</b>	<b>Sasharomyces Cerevisiae</b>
0.1µm	>7/cm <sup>2</sup>	N/A	N/A
0.2µm	>7/cm <sup>2</sup>	N/A	N/A
0.45µm	N/A	>7/cm <sup>2</sup>	>7/cm <sup>2</sup>
0.65µm	N/A	>4/cm <sup>2</sup>	>7/cm <sup>2</sup>
1.2µm	N/A	N/A	>7/cm <sup>2</sup>

Log Reduction Values are calculated using the following formula:  $LRV = \log_{10} \left( \frac{\text{total number of organisms entering the filter}}{\text{total number of organisms exiting the filter}} \right)$

## Ordering Information

<b>APSBR</b>	<b>Removal</b>	<b>End Cap</b>	<b>Nominal Length</b>	<b>Seal Material</b>	<b>-F</b>
[Φ69]	<b>0010</b> =0.1 µm	<b>DOE</b> =Double open end	<b>05</b> = 5"	<b>S</b> =Silicone	
	<b>0022</b> =0.22 µm	<b>HTC</b> =222 o-ring/flat (PBT Insert)	<b>10</b> =10"	<b>E</b> =EPDM	
	<b>0045</b> =0.45 µm	<b>HTF</b> =222 o-ring/fin (PBT Insert)	<b>20</b> =20"	<b>V</b> =Viton	
	<b>0065</b> =0.65 µm	<b>HSF</b> =226 o-ring/fin (PBT Insert)	<b>30</b> =30"		
	<b>0080</b> =0.8 µm	<b>SSF</b> =226 o-ring/fin (SS insert)	<b>40</b> =40"		
	<b>0120</b> =1.2 µm	<b>SSC</b> =226 o-ring/flat (SS insert)			
		<b>STF</b> =222 o-ring/fin (SS insert, 3 Tabs)			

# Hall Pyke



[www.hallpyke.ie](http://www.hallpyke.ie)

Head Office, 3A Sunbury Industrial Estate, Ballymount Road, Walkinstown, Dublin 12.  
T. +353-1-4501411 F. +353-1-4507960 E. [info@hallpyke.ie](mailto:info@hallpyke.ie) W. [www.hallpyke.ie](http://www.hallpyke.ie)